



INSPEC International Limited
56 Leslie Hough Way
Salford, Greater Manchester
M6 6AJ, United Kingdom

Tel: +44 (0) 16 17 37 06 99
Fax: +44 (0) 16 17 36 01 01
Email: equipment@inspec-international.com

Acoustic test fixture

Description:

As fully described in ISO 24869-3 13819-2, the acoustic test fixture is designed to accommodate a 1" microphone and pre-amplifier. It is provided with an isolation cup and (optionally) an adaptor to accommodate earmuffs mounted on an industrial safety helmet.

Services required: None

Approximate packed size & weight: 40 x 30 x 30cm : 8kg

Relevant standard: ISO 4869-3, para 5.1 and fig 1.

Delivery: Approximately 6 weeks from receipt of order and deposit.

Effective: May 2008

