

Stolle Apertometer

Description:

The apertometer comprises two major assemblies:

- A headform with light bulbs located at the eye positions
- A graduated, translucent, thermo-formed plastic bowl of nominal 610mm diameter.

These are fixed to a suitable base plate to enable bench mounting.

All as described in EN 136:1998, clause 8.17.

Not included, but recommended, is a planimeter, suitable for measuring the field of vision chart produced using the apertometer above. Please refer to separate specification.

Services required:

Bench mounted. 230 or 115 vac, 50 / 60Hz mains electricity. Please specify.

Approximate packed size & weight: 90 x 60 x 50 cm : 25 kg

Relevant standards: EN 136:1998, clause 8.17

Delivery: Approximately 12 weeks from receipt of order and deposit.

Effective: May 2008

