



INSPEC International Limited
56 Leslie Hough Way
Salford, Greater Manchester
M6 6AJ, United Kingdom

Tel: +44 (0) 16 17 37 06 99
Fax: +44 (0) 16 17 36 01 01
Email: equipment@inspec-international.com

Dolomite dust test equipment

Description:

The equipment comprises the chamber, as described in EN143, together with an automatic dust dispenser, air mover to draw air through the chamber, an electronic dust deposition indicator, digital manometer and chart recorder. The sample is supported in the dust stream by means of a suitable fixture.

Services required:

Floor mounted. 230/115vac, 50/60Hz mains electricity. Please specify.
Instrument quality compressed air.

Approximate packed size & weight: To be advised.

Relevant standards: ENs 143 / 149 etc

Delivery: Approximately 6 to 8 weeks from receipt of order and deposit.

Effective: May 2008

