

## **Drop test plate**

### ***Description:***

A mild steel plate, 500 x 500 x 12 mm, edged & chamfered, flash ground top and bottom and zinc plated.

### ***Services required:***

None. To be placed on a firm, level floor.

### ***Approximate packed size & weight:***

70 x 70 x 30 cm : 7kg

### ***Relevant standard:***

BS EN 13819-1:2002 (and EN 352 Parts 1, 2 and 3)

***Delivery:*** 3 weeks from receipt of order and deposit.

Effective: May 2008

