

Electrical Insulation Test Equipment

Description:

The three tests in the EN 397 and EN 443 standards are covered.

Two tests utilise the standard aluminium headform and the other uses a special water tank with electrical probes built in.

A self-contained electrical tester is then used to measure the leakage current at specified voltages and at various points on the helmet in different configurations.

Services required:

Bench mounted. Either 230 or 115vac, 50/60Hz mains electricity. Please specify.

Approximate packed size & weight:

50 x 80 x 40 cm : 25 kg

Relevant standard:

EN 397 : 1995, paras 6.10.1, 6.10.2 and 6.10.3

EN443 : 2008, paras 5.14.1, 5.14.2 and 5.14.3

EN 13087-8: 2000

Delivery: Approximately 8 weeks from receipt of order and deposit.

Effective: May 2008

