

## **Set of 3 test headforms for helmet-mounted earmuffs**

### ***Description:***

Each of the three headforms (size codes B, J and N) comprises a cast epoxy resin “half-head”, secured between two parallel, transparent acrylic plates. On each plate, at the points corresponding to the centre of the ear openings, is located a threaded insert to enable the attachment of an artificial pinna at the appropriate point. The headform and plates assembly is attached to a hardwood baseplate.

One pair of pinnae and attachment screws are included.

### ***Services required:***

Bench mounted. No other services required.

### ***Approximate packed size & weight:***

80 x 35 x 40 cm : 10 kg

### ***Relevant standard:***

EN13819-1:2001, para 4.2.2.2 and fig 4

***Delivery:*** Approximately 6 weeks from receipt of order and deposit.

Effective: May 2008

