



INSPEC International Limited
56 Leslie Hough Way
Salford, Greater Manchester
M6 6AJ, United Kingdom

Tel: +44 (0) 16 17 37 06 99
Fax: +44 (0) 16 17 36 01 01
Email: equipment@inspec-international.com

Acoustic instrumentation for insertion loss test

Description:

The acoustic instrumentation package includes a Bruel & Kjaer 1" microphone and pre-amplifier (for the ATF), Pulse system, power amplifier (for the signal generation and analysis) and microphone calibrator - fully described in, and complying with, ISO 4869-3 & EN 13819-2.

All necessary connectors and cables are included.

A laptop computer is pre-loaded with the necessary software to run the insertion loss tests and present the results in a pre-formatted test report. Data may be saved to an Excel spreadsheet.

Services required:

Bench mounted. Either 230 or 115vac, 50/60Hz mains electricity. Please specify.

Approximate packed size & weight: To be advised

Relevant standards: ISO 4869-3 and EN 13819-2

Delivery: Approximately 6 weeks from receipt of order and deposit.

Effective: May 2008

