

Ignitability Test Equipment

Description:

Comprises a 300mm long by 6mm diameter steel rod, with a thermocouple inside, whose junction is located 20mm from the end of the rod to be heated. The thermocouple is connected to an electronic digital thermometer.

A hand held propane burner is suggested for heating the end of the rod.

Services required:

Battery power for thermometer. It is suggested that the burner should be purchased locally, to ensure a continuing supply of replacement gas cartridges.

Approximate packed size & weight:

40 x 15 x 12 cm : 2 kg

Relevant standard:

EN 168 : 2001, para 7.1. EN 13819-1: 2002, para 4.13.2

Delivery: 4 weeks from receipt of order and deposit.

Effective: May 2008

