



INSPEC International Limited  
56 Leslie Hough Way  
Salford, Greater Manchester  
M6 6AJ, United Kingdom

Tel: +44 (0) 16 17 37 06 99  
Fax: +44 (0) 16 17 36 01 01  
Email: [equipment@inspec-international.com](mailto:equipment@inspec-international.com)

## Polarizing Axis Test Equipment

### ***Description:***

Equipment exactly as specified in 6.1.3.2.1 of EN1836 : 2001 and Annex D.1 of EN 172 : 2000

### ***Services required:***

Bench mounted. 230vac, 50Hz mains electricity

### ***Approximate packed size & weight:***

40cm x 40cm x 40cm : 5kg

### ***Relevant standards:***

EN 172 : 2000  
EN1836 : 2001

***Delivery:*** 4 weeks from receipt of order and deposit.

Effective: May 2008

