

## **Falling Sand Test Equipment (Resistance to Surface Damage by Fine Particles)**

### ***Description:***

The falling sand apparatus, as specified in EN 168:2002, para 15.1 and figures 8 & 9.

### ***Services required:***

Floor mounted. 230/115vac, 50/60Hz mains electricity. Please specify.

Consumables include 3 kg of specified sieved sand per 10 tests, and samples in reference glass and PMMA. The quotation includes sufficient of these consumables for 50 tests.

### ***Approximate packed size & weight:***

250 x 50 x 50 cm : 65 kg

### ***Relevant standard:***

BS EN 168:2002, para 15.1

***Delivery:*** 12 weeks from receipt of order and deposit.

Effective: May 2008

